

# Laser Esthetics

SmoothLase = Wrinkle Reduction  
LipLase = Natural Lip Enhancement  
Laser Skin Resurfacing



Start Looking & Feeling Years Younger Naturally

- See Instant Results
- No Injections
- No Recovery Time
- Improve Skin's Color & Tone
- Long Lasting Results You Can See & Feel
- Stimulates Your Body's Own Production of Collagen
- Helps You Look & Feel Better Naturally

twissDENTAL

## Why Twiss Dental is your best option for SmoothLase™, LipLase™, & Skin Resurfacing

Our health care providers have taken the advanced training necessary to perform these intricate procedures. We can help restore your youthfulness by stimulating your body's own natural collagen production . Call now to see if we can help.

### Twiss Dental

Tyler Twiss, DMD | MBA

4185 E Wildcat Reserve Pkwy #110  
Highlands Ranch, CO 80126

720.770.4990

[twissdental.com](http://twissdental.com)

 [twissdentalco](https://www.facebook.com/twissdentalco)

 [twissdental](https://www.instagram.com/twissdental)

**Fotona**  
Ultra Performance Lasers™

SmoothLase - Before & After

LipLase - Before & After


Laser Skin Resurfacing - Before & After

"I love the smooth, natural look of younger skin that SmoothLase gave me."

Denise C.

photo courtesy of Franco Paciolla, MD

Treating the entire oral cavity can result in the lifting & tightening of the nasolabial folds, lip lines, jowls, & even a saggy chin.



SmoothLase = Wrinkle Reduction (AKA "Inside-Out Facelift")

## How SmoothLase & LipLase Work

SmoothLase and LipLase use the LightWalker™ laser to generate the formation of new and more elastic collagen. The "Smooth Mode" pulse creates a non-ablative heat which initiates new collagen formation. The LightWalker™ handpiece is inserted intraorally with no direct contact to the inside of your mouth. The laser emits pulses of light to the inside of your mouth without damage to your tissues. In LipLase, the laser delivers targeted heat pulses into the mucosa of your lips to stimulate new collagen formation resulting in fuller, more defined and enhanced lips. A consultation will determine if you need more than one treatment with additional treatments spaced over a few months.

The Erbium YAG laser transforms old, saggy Category III collagen back to firm and tighter Category I collagen. This naturally lifts and tightens the tissue, reversing the signs of aging in most cases. Treatment is comprised of 3-4 painless visits. The Erbium YAG laser also stimulates new collagen to form, preserving the results.



Category III Collagen Stimulation



Laser Stimulated Change



Category I Collagen



Erbium Neocollagenesis

